



PPM1F Antibody, HRP conjugated

Product Code	CSB-PA018493LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49593
Immunogen	Recombinant Human Protein phosphatase 1F protein (1-76AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein phosphatase 1F (EC 3.1.3.16) (Ca(2+)/calmodulin-dependent protein kinase phosphatase) (CaM-kinase phosphatase) (CaMKPase) (Partner of PIX 2) (Protein fem-2 homolog) (hFem-2), PPM1F, KIAA0015 POPX2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PPM1F
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.