



PTPN12 Antibody, HRP conjugated

Product Code	CSB-PA019026LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05209
Immunogen	Recombinant Human Tyrosine-protein phosphatase non-receptor type 12 protein (539-714AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein phosphatase non-receptor type 12 (EC 3.1.3.48) (PTP-PEST) (Protein-tyrosine phosphatase G1) (PTPG1), PTPN12
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PTPN12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.