



# RECK Antibody, Biotin conjugated

Product Code	CSB-PA019536LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95980
Immunogen	Recombinant Human Reversion-inducing cysteine-rich protein with Kazal motifs protein (23-212AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Reversion-inducing cysteine-rich protein with Kazal motifs (hRECK) (Suppressor of tumorigenicity 15 protein), RECK, ST15
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	RECK