



RECK Antibody, Biotin conjugated

Product Code	CSB-PA019536LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95980
Immunogen	Recombinant Human Reversion-inducing cysteine-rich protein with Kazal motifs protein (23-212AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Reversion-inducing cysteine-rich protein with Kazal motifs (hRECK) (Suppressor of tumorigenicity 15 protein), RECK, ST15
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	RECK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.