



SCARA5 Antibody, FITC conjugated

Product Code	CSB-PA020766LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZMJ2
Immunogen	Recombinant Human Scavenger receptor class A member 5 protein (303-428AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Scavenger receptor class A member 5 (Scavenger receptor hlg), SCARA5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SCARA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.