



SCNN1B Antibody, HRP conjugated

Product Code	CSB-PA020849LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51168
Immunogen	Recombinant Human Amiloride-sensitive sodium channel subunit beta protein (554-640AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Amiloride-sensitive sodium channel subunit beta (Beta-NaCH) (Epithelial Na(+) channel subunit beta) (Beta-ENaC) (ENaCB) (Nonvoltage-gated sodium channel 1 subunit beta) (SCNEB), SCNN1B
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SCNN1B