



SLC11A1 Antibody, HRP conjugated

Product Code	CSB-PA021380LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49279
Immunogen	Recombinant Human Natural resistance-associated macrophage protein 1 protein (1-58AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Natural resistance-associated macrophage protein 1 (NRAMP 1) (Solute carrier family 11 member 1), SLC11A1, LSH NRAMP NRAMP1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SLC11A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.