



# SLC11A1 Antibody, Biotin conjugated

|                     |   |
|---------------------|---|
| Product Code        | CSB-PA021380LD01HU  |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| Uniprot No.         | P49279  |
| Immunogen           | Recombinant Human Natural resistance-associated macrophage protein 1 protein (1-58AA)                                       |
| Raised In           | Rabbit  |
| Species Reactivity  | Human   |
| Tested Applications | ELISA   |
| Form                | Liquid  |
| Conjugate           | Biotin  |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4  |
| Purification Method | >95%, Protein G purified  |
| Isotype             | IgG   |
| Clonality           | Polyclonal  |
| Alias               | Natural resistance-associated macrophage protein 1 (NRAMP 1) (Solute carrier family 11 member 1), SLC11A1, LSH NRAMP NRAMP1 |
| Immunogen Species   | Homo sapiens (Human)  |
| Research Area       | Immunology  |
| Target Names        | SLC11A1   |