



SLC12A2 Antibody, HRP conjugated

Product Code	CSB-PA021383LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55011
Immunogen	Recombinant Human Solute carrier family 12 member 2 protein (110-282AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 12 member 2 (Basolateral Na-K-Cl symporter) (Bumetanide-sensitive sodium-(potassium)-chloride cotransporter 1), SLC12A2, NKCC1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SLC12A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.