



SNAP23 Antibody, FITC conjugated

Product Code	CSB-PA021871LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00161
Immunogen	Recombinant Human Synaptosomal-associated protein 23 protein (90-196AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Synaptosomal-associated protein 23 (SNAP-23) (Vesicle-membrane fusion protein SNAP-23), SNAP23
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SNAP23