



SREBF1 Antibody, Biotin conjugated

Product Code	CSB-PA022657LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36956
Immunogen	Recombinant Human Sterol regulatory element-binding protein 1 protein (38-201AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sterol regulatory element-binding protein 1 (SREBP-1) (Class D basic helix-loop-helix protein 1) (bHLHd1) (Sterol regulatory element-binding transcription factor 1) [Cleaved into: Processed sterol regulatory element-binding protein 1], SREBF1, BHLHD1 SREBP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SREBF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.