



TAAR5 Antibody, Biotin conjugated

Product Code	CSB-PA023057LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14804
Immunogen	Recombinant Human Trace amine-associated receptor 5 protein (1-34AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trace amine-associated receptor 5 (TaR-5) (Trace amine receptor 5) (hTaar5) (Putative neurotransmitter receptor), TAAR5, PNR
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TAAR5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.