



TCP1 Antibody, Biotin conjugated

Product Code	CSB-PA023320LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17987
Immunogen	Recombinant Human T-complex protein 1 subunit alpha protein (406-506AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-complex protein 1 subunit alpha (TCP-1-alpha) (CCT-alpha), TCP1, CCT1 CCTA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TCP1