



TNFAIP2 Antibody, FITC conjugated

Product Code	CSB-PA023957LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q03169
Immunogen	Recombinant Human Tumor necrosis factor alpha-induced protein 2 protein (1-93AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor alpha-induced protein 2 (TNF alpha-induced protein 2) (Primary response gene B94 protein), TNFAIP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TNFAIP2