



TPP1 Antibody, Biotin conjugated

Product Code	CSB-PA024113LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14773
Immunogen	Recombinant Human Tripeptidyl-peptidase 1 protein (371-480AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripeptidyl-peptidase 1 (TPP-1) (EC 3.4.14.9) (Cell growth-inhibiting gene 1 protein) (Lysosomal pepstatin-insensitive protease) (LPIC) (Tripeptidyl aminopeptidase) (Tripeptidyl-peptidase I) (TPP-I), TPP1, CLN2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TPP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.