



ZFP36L2 Antibody, HRP conjugated

Product Code	CSB-PA026450LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47974
Immunogen	Recombinant Human mRNA decay activator protein ZFP36L2 protein (219-321AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	mRNA decay activator protein ZFP36L2 (Butyrate response factor 2) (EGF-response factor 2) (ERF-2) (TPA-induced sequence 11d) (Zinc finger protein 36, C3H1 type-like 2) (ZFP36-like 2), ZFP36L2, BRF2 ERF2 RNF162C TIS11D
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ZFP36L2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.