



PNMT Antibody, FITC conjugated

Product Code	CSB-PA02939C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11086
Immunogen	Recombinant Human Phenylethanolamine N-methyltransferase protein (1-176AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phenylethanolamine N-methyltransferase (PNMTase) (EC 2.1.1.28) (Noradrenaline N-methyltransferase), PNMT, PENT
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PNMT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.