



HSPB1 Antibody, Biotin conjugated

Product Code	CSB-PA10869D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04792
Immunogen	Recombinant Human Heat shock protein beta-1 protein (108-205AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock protein beta-1 (HspB1) (28 kDa heat shock protein) (Estrogen-regulated 24 kDa protein) (Heat shock 27 kDa protein) (HSP 27) (Stress-responsive protein 27) (SRP27), HSPB1, HSP27 HSP28
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSPB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.