



CCKBR Antibody, HRP conjugated

Product Code	CSB-PA12619B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32239
Immunogen	Recombinant Human Gastrin/cholecystokinin type B receptor protein (243-333AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gastrin/cholecystokinin type B receptor (CCK-B receptor) (CCK-BR) (Cholecystokinin-2 receptor) (CCK2-R), CCKBR, CCKRB
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CCKBR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.