



IL10RB Antibody, HRP conjugated

Product Code	CSB-PA600091LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08334
Immunogen	Recombinant Human Interleukin-10 receptor subunit beta protein (79-219AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-10 receptor subunit beta (IL-10 receptor subunit beta) (IL-10R subunit beta) (IL-10RB) (Cytokine receptor class-II member 4) (Cytokine receptor family 2 member 4) (CRF2-4) (Interleukin-10 receptor subunit 2) (IL-10R subunit 2) (IL-10R2) (CD antigen CDw210b), IL10RB, CRFB4 D21S58 D21S66
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL10RB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.