



MTMR2 Antibody, Biotin conjugated

Product Code	CSB-PA613418LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13614
Immunogen	Recombinant Human Myotubularin-related protein 2 protein (1-173AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myotubularin-related protein 2 (Phosphatidylinositol-3,5-bisphosphate 3-phosphatase) (EC 3.1.3.95) (Phosphatidylinositol-3-phosphate phosphatase) (EC 3.1.3.64), MTMR2, KIAA1073
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MTMR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.