



SEMA3A Antibody, FITC conjugated

Product Code	CSB-PA613506LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14563
Immunogen	Recombinant Human Semaphorin-3A protein (357-445AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Semaphorin-3A (Semaphorin III) (Sema III), SEMA3A, SEMAD
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SEMA3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.