



LAMB3 Antibody, HRP conjugated

Product Code	CSB-PA615717LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13751
Immunogen	Recombinant Human Laminin subunit beta-3 protein (256-397AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Laminin subunit beta-3 (Epiligrin subunit bata) (Kalinin B1 chain) (Kalinin subunit beta) (Laminin B1k chain) (Laminin-5 subunit beta) (Nicein subunit beta), LAMB3, LAMNB1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LAMB3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.