



GRM2 Antibody, FITC conjugated

Product Code	CSB-PA621775LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14416
Immunogen	Recombinant Human Metabotropic glutamate receptor 2 protein (93-221AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metabotropic glutamate receptor 2 (mGluR2), GRM2, GPRC1B MGLUR2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.