



RAB3GAP1 Antibody, Biotin conjugated

Product Code	CSB-PA621861LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15042
Immunogen	Recombinant Human Rab3 GTPase-activating protein catalytic subunit protein (527-653AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rab3 GTPase-activating protein catalytic subunit (RAB3 GTPase-activating protein 130 kDa subunit) (Rab3-GAP p130) (Rab3-GAP), RAB3GAP1, KIAA0066 RAB3GAP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RAB3GAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.