



IQCB1 Antibody, HRP conjugated

Product Code	CSB-PA621862LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15051
Immunogen	Recombinant Human IQ calmodulin-binding motif-containing protein 1 protein (445-598AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	IQ calmodulin-binding motif-containing protein 1 (Nephrocystin-5) (p53 and DNA damage-regulated IQ motif protein) (PIQ), IQCB1, KIAA0036 NPHP5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IQCB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.