



IQCB1 Antibody, FITC conjugated

Product Code	CSB-PA621862LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15051
Immunogen	Recombinant Human IQ calmodulin-binding motif-containing protein 1 protein (445-598AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	IQ calmodulin-binding motif-containing protein 1 (Nephrocystin-5) (p53 and DNA damage-regulated IQ motif protein) (PIQ), IQCB1, KIAA0036 NPHP5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IQCB1