



DUSP6 Antibody, Biotin conjugated

Product Code	CSB-PA621973LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16828
Immunogen	Recombinant Human Dual specificity protein phosphatase 6 protein (100-206AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity protein phosphatase 6 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity protein phosphatase PYST1) (Mitogen-activated protein kinase phosphatase 3) (MAP kinase phosphatase 3) (MKP-3), DUSP6, MKP3 PYST1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DUSP6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.