



UBAP2L Antibody, HRP conjugated

Product Code	CSB-PA622754LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14157
Immunogen	Recombinant Human Ubiquitin-associated protein 2-like protein (104-252AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-associated protein 2-like (Protein NICE-4), UBAP2L, KIAA0144 NICE4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UBAP2L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.