



PCBP2 Antibody, Biotin conjugated

Product Code	CSB-PA622995LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15366
Immunogen	Recombinant Human Poly(rC)-binding protein 2 protein (171-293AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Poly(rC)-binding protein 2 (Alpha-CP2) (Heterogeneous nuclear ribonucleoprotein E2) (hnRNP E2), PCBP2
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	PCBP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.