



PRSS8 Antibody, Biotin conjugated

Product Code	CSB-PA623091LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16651
Immunogen	Recombinant Human Prostasin protein (185-319AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prostasin (EC 3.4.21.-) (Channel-activating protease 1) (CAP1) (Serine protease 8) [Cleaved into: Prostasin light chain; Prostasin heavy chain], PRSS8
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PRSS8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.