



LAMA3 Antibody, Biotin conjugated

Product Code	CSB-PA623096LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16787
Immunogen	Recombinant Human Laminin subunit alpha-3 protein (1570-1796AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Laminin subunit alpha-3 (Epiligrin 170 kDa subunit) (E170) (Epiligrin subunit alpha) (Kalinin subunit alpha) (Laminin-5 subunit alpha) (Laminin-6 subunit alpha) (Laminin-7 subunit alpha) (Nicein subunit alpha), LAMA3, LAMNA
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LAMA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.