



TRDMT1 Antibody, HRP conjugated

Product Code	CSB-PA675694LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14717
Immunogen	Recombinant Human tRNA (cytosine(38)-C(5))-methyltransferase protein (170-269AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	tRNA (cytosine(38)-C(5))-methyltransferase (EC 2.1.1.204) (DNA (cytosine-5)-methyltransferase-like protein 2) (Dnmt2) (DNA methyltransferase homolog HsaIIP) (DNA MTase homolog HsaIIP) (M.HsaIIP) (PuMet), TRDMT1, DNMT2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TRDMT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.