



TRDMT1 Antibody, HRP conjugated

Product Code	CSB-PA675694LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14717
Immunogen	Recombinant Human tRNA (cytosine(38)-C(5))-methyltransferase protein (170-269AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	tRNA (cytosine(38)-C(5))-methyltransferase (EC 2.1.1.204) (DNA (cytosine-5)-methyltransferase-like protein 2) (Dnmt2) (DNA methyltransferase homolog HsaIIP) (DNA MTase homolog HsaIIP) (M.HsaIIP) (PuMet), TRDMT1, DNMT2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TRDMT1