



MYBPC1 Antibody, Biotin conjugated

Product Code	CSB-PA717389LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00872
Immunogen	Recombinant Human Myosin-binding protein C, slow-type protein (618-779AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin-binding protein C, slow-type (Slow MyBP-C) (C-protein, skeletal muscle slow isoform), MYBPC1, MYBPCS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYBPC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.