



MYBPC1 Antibody, Biotin conjugated

Product Code	CSB-PA717389LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00872
Immunogen	Recombinant Human Myosin-binding protein C, slow-type protein (618-779AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin-binding protein C, slow-type (Slow MyBP-C) (C-protein, skeletal muscle slow isoform), MYBPC1, MYBPCS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYBPC1