



GPRC6A Antibody, Biotin conjugated

Product Code	CSB-PA722571LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T6X5
Immunogen	Recombinant Human G-protein coupled receptor family C group 6 member A protein (403-581AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	G-protein coupled receptor family C group 6 member A (hGPRC6A) (G-protein coupled receptor GPCR33) (hGPCR33), GPRC6A
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GPRC6A