



SLC22A25 Antibody, Biotin conjugated

Product Code	CSB-PA740875LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6T423
Immunogen	Recombinant Human Solute carrier family 22 member 25 protein (31-145AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 22 member 25 (Organic anion transporter UST6), SLC22A25, UST6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC22A25
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.