



TRAF7 Antibody, Biotin conjugated

Product Code	CSB-PA757671LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6Q0C0
Immunogen	Recombinant Human E3 ubiquitin-protein ligase TRAF7 protein (1-90AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase TRAF7 (EC 2.3.2.27) (RING finger and WD repeat-containing protein 1) (RING finger protein 119) (RING-type E3 ubiquitin transferase TRAF7) (TNF receptor-associated factor 7), TRAF7, RFWD1 RNF119
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRAF7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.