



GDF7 Antibody, FITC conjugated

Product Code	CSB-PA768769LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z4P5
Immunogen	Recombinant Human Growth/differentiation factor 7 protein (322-450AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth/differentiation factor 7 (GDF-7), GDF7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GDF7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.