



DPP9 Antibody, FITC conjugated

Product Code	CSB-PA769770LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86TI2
Immunogen	Recombinant Human Dipeptidyl peptidase 9 protein (289-437AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dipeptidyl peptidase 9 (DP9) (EC 3.4.14.5) (Dipeptidyl peptidase IV-related protein 2) (DPRP-2) (Dipeptidyl peptidase IX) (DPP IX) (Dipeptidyl peptidase-like protein 9) (DPLP9), DPP9, DPRP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DPP9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.