



USP48 Antibody, Biotin conjugated

Product Code	CSB-PA774808LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UV5
Immunogen	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 48 protein (358-508AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase 48 (EC 3.4.19.12) (Deubiquitinating enzyme 48) (Ubiquitin thioesterase 48) (Ubiquitin-specific peptidase 48) (Ubiquitin-specific protease 48) (Ubiquitin-specific-processing protease 48), USP48, USP31
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	USP48
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.