



PNPLA6 Antibody, Biotin conjugated

Product Code	CSB-PA809011LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IY17
Immunogen	Recombinant Human Neuropathy target esterase protein (345-493AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuropathy target esterase (EC 3.1.1.5) (Patatin-like phospholipase domain-containing protein 6), PNPLA6, NTE
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PNPLA6