



NPAS4 Antibody, HRP conjugated

Product Code	CSB-PA814204LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IUM7
Immunogen	Recombinant Human Neuronal PAS domain-containing protein 4 protein (445-585AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuronal PAS domain-containing protein 4 (Neuronal PAS4) (Class E basic helix-loop-helix protein 79) (bHLHe79) (HLH-PAS transcription factor NXF) (PAS domain-containing protein 10), NPAS4, BHLHE79 NXF PASD10
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NPAS4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.