



IRF2BP1 Antibody, HRP conjugated

Product Code	CSB-PA815541LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IU81
Immunogen	Recombinant Human Interferon regulatory factor 2-binding protein 1 protein (340-499AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon regulatory factor 2-binding protein 1 (IRF-2-binding protein 1) (IRF-2BP1) (Probable E3 ubiquitin-protein ligase IRF2BP1) (EC 2.3.2.27) (Probable RING-type E3 ubiquitin transferase IRF2BP1), IRF2BP1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IRF2BP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.