



GPBAR1 Antibody, HRP conjugated

Product Code	CSB-PA819471LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDU6
Immunogen	Recombinant Human G-protein coupled bile acid receptor 1 protein (283-330AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	G-protein coupled bile acid receptor 1 (G-protein coupled receptor GPCR19) (hGPCR19) (Membrane-type receptor for bile acids) (M-BAR) (hBG37) (BG37), GPBAR1, TGR5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPBAR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.