



CAMKK2 Antibody, Biotin conjugated

Product Code	CSB-PA822299LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RR4
Immunogen	Recombinant Human Calcium/calmodulin-dependent protein kinase kinase 2 protein (5-148AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium/calmodulin-dependent protein kinase kinase 2 (CaM-KK 2) (CaM-kinase kinase 2) (CaMKK 2) (EC 2.7.11.17) (Calcium/calmodulin-dependent protein kinase kinase beta) (CaM-KK beta) (CaM-kinase kinase beta) (CaMKK beta), CAMKK2, CAMKKB KIAA0787
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CAMKK2