



CAMKK2 Antibody, Biotin conjugated

Product Code	CSB-PA822299LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RR4
Immunogen	Recombinant Human Calcium/calmodulin-dependent protein kinase kinase 2 protein (5-148AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium/calmodulin-dependent protein kinase kinase 2 (CaM-KK 2) (CaM-kinase kinase 2) (CaMKK 2) (EC 2.7.11.17) (Calcium/calmodulin-dependent protein kinase kinase beta) (CaM-KK beta) (CaM-kinase kinase beta) (CaMKK beta), CAMKK2, CAMKKB KIAA0787
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CAMKK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.