



LGI2 Antibody, Biotin conjugated

Product Code	CSB-PA822674LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N0V4
Immunogen	Recombinant Human Leucine-rich repeat LGI family member 2 protein (194-310AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat LGI family member 2 (LGI1-like protein 2) (Leucine-rich glioma-inactivated protein 2), LGI2, KIAA1916 LGIL2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LGI2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.