



NUP93 Antibody, FITC conjugated

Product Code	CSB-PA822683LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N1F7
Immunogen	Recombinant Human Nuclear pore complex protein Nup93 protein (14-180AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear pore complex protein Nup93 (93 kDa nucleoporin) (Nucleoporin Nup93), NUP93, KIAA0095
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NUP93